



Wataru Obuchi
Laboratory Manager
Sydney, Australia

“Creating the SomnoDent Fit”

Practice Details:

NOTE: SomnoMed recommends that all patients are diagnosed by a registered Sleep Physician before undergoing treatment.

SomnoDent®

SLEEP APNEA APPLIANCE



Why SomnoDent?

Permits normal mouth opening and closing

Eliminates the feeling of claustrophobia; allows speech and drinking during night.

Fits comfortably in mouth

Discreet when worn, the appliance provides normal tongue space and permits full lip seal.

Incrementally adjustable

Fine tune protrusive adjustments on both sides independently for optimum positioning; provides more comfort and improves patient tolerance.

Stable, flat occlusion

Even distribution of forces over the two plates, providing better protection for teeth.

Clinically demonstrated

Demonstrated to work in most cases of mild to moderate OSA and snoring.

Quality construction

ISO certification, CE registered and Swissmedic, FDA cleared.

Customisable to suit patients' needs

Devices are available with additional features which may suit bruxers and mouth breathers.

Patented and registered

Europe, Switzerland, Canada, United States, Australia, South Korea, Japan and Hong Kong

SomnoMed Limited

Level 3, 20 Clarke Street
Crows Nest NSW 2065
phone: 1800 445 660
fax +61 2 9467 0467

www.somnomed.com.au

Product Guide

- *Patented*
- *Clinically Proven*
- *Patient Compliant*



SomnoDent®
FLEX



Made from our unique SMH BFlex material and acrylic that provides superior retention and comfort. SMH BFlex is a soft polymer that is heat-cured, packed and pressed utilising an exclusive chemical formula that facilitates soft to hard polymer bonding at the layers interface. The result is a fit that is unmatched and unlike anything seen before in a resilient material, securely anchoring appliance with a snug fit – even on short teeth with little undercut.

With an expected life of 3-5 years, it is suitable for all dentitions. SMH BFlex does not discolour or attract odours. Patients and dentists may prefer this model, particularly where crown and bridge work is present.

SomnoDent®
CLASSIC



The original SomnoDent, made from the highest quality acrylic that will not discolour, the Classic has been specifically designed to offer solidity, comfort and quality. The Classic is retained using ball clasps and has a life expectancy of 3-5 years.

SomnoDent®
EDENT



Designed for either upper edentulous or completely edentulous patients, the Edent and Edent Implant SomnoDents can be made from either SMH BFlex or acrylic materials.

The Edent upper splint covers the palate in the same way as a full upper denture secured by suction and made to oppose a lower splint. Due to the unique two piece coupling mechanism there is no downward pull on the upper splint. This advantage allows edentulous patients to enjoy the benefits of a SomnoDent®.

For completely edentulous patients, SomnoMed recommends the Edent implant.

 Gold Standard



All SomnoDents now come with a 2 year warranty
*see website for terms and conditions

SomnoDent®
SLEEP APNEA APPLIANCE

- Safe and effective
- Clinically validated
- Custom-Fitted
- Permits normal mouth opening
- Allows speech and drinking
- Warranty against breakages
- Full lip-seal
- Easily Adjustable